Cryotherapy

Targeted surface cooling of tissue



TECHNOLOGY FOR THERAPY



FRIGOSTREAM

Environmentally friendly cryotherapy in a compact design

The compact FRIGOSTREAM cryotherapy unit offers optimised performance in a minimum of space and is available in five different versions. It generates cold air with an adjustable current of air of up to 1500 l/min and thus enables targeted cooling of body areas with the recognised advantages over CO₂ and ice packs. If used correctly the risk of cold damage is excluded.

FRIGOSTREAM operates with an energy source that costs nothing - air - and is therefore always ready for use and has no resource replenishing problems. Like an air conditioning system, the ambient air is cooled down and blown onto the skin through a treatment tube without direct contact at temperatures of up to -60 °C (see technical data). Various tubes and nozzles permit easy targeted treatment.

The de-icing function, which requires no water emptying on the user's part (with the exception of C200/ C600), is particularly user-friendly. The FRIGOSTREAM line refrigerant loop is a closed system and requires no maintenance whatsoever.



Pain therapy with cooled air



Cryotherapy

Cryotherapy is based on surface cooling and thus local reduction of tissue temperature. Cold air therapy uses a continuously cooled stream of air within a constant temperature range, which results in a series of thermoregulatory effects that differ from other cold application forms (e.g. ice or cold packs) and are desirable for various indications, including acute and chronic pain (especially in sport physiotherapy) and rheumatism. Cold air therapy is also used to cool the skin surface during laser treatment.

In detail, there are the following specific physiological effects of cryotherapy with cold air flow:

- » Reduction in the sensitivity of sensory nerves
- » Reduction in nerve conduction
- » Immediate vasoconstrictive effect

SPECIAL FEATURES

Cryotherapie

	Cold air from –32 °C to –60 °C			
	Adjustable stream of air up to 1500 l/min.			
	Very short cooling down phase			
	De-icing function Maintenance-free refrigerant loop			
	Standard slotted nozzle			
(1)	Nozzle attachments in three sizes for effective treatment			
2	Flexible, double-wall tube (1.70 m long)			
3	Therapy arm for fast and easy placement of the tube over a clip-on system (option)			
	FRIGOSTREAM is also available with a special smaller dia-			

meter tube to maintain cooling at a constant level over an

extended period for thermal laser treatments

GENERAL FEATURES

Intuitive operation Touch screen







TECHNICAL DATA

FRIGOSTREAM, FRIGOSTREAM Premium and Turbo:

Protection class	1, type B			
Power connection	230 V ±10 % or 110 V ±10 %*			
Mains frequency	50 – 60 Hz			
Power consumption	Ø 650 W (500 W**, 850 W***)			
Cooling capacity	up to -32 °C (-40 °C***)			
Air capacity	350 — 1500 l/min. (350 — 1200 l/min**)			
Dimensions (W x H x D)	360 x 930 x 480 mm			
Weight	53 kg			
*Turbo only 230 V ±10 %	**Premium ***Turbo			

FRIGOSTREAM C200 und C600:

Protection class	1, type B	
Power connection	230 V ±10 %	
Mains frequency	50 – 60 Hz	
Power consumption	900 W (1400 W*)	
Cooling capacity	up to -40 °C (-60 °C*)	
Air capacity	350 – 1500 l/min.	
Dimensions (W x H x D)	380 x 1020 x 540 mm	
Weight	65 kg (80 kg*)	
	*C 600	

STANDARD ACCESSORIES

[1]	Nozzle attachment 5 mm
[1]	Nozzle attachment 10 mm
[1]	Nozzle attachment 15 mm
[1]	Operating instructions
[1]	Slotted nozzle
[1]	Tube

ADDRESS

PHONE

E-MAIL

PHYSIOMED ELEKTROMEDIZIN AG	+49 (0) 91 26 / 25 87-0	info@physiomed.de
Hutweide 10		
91220 Schnaittach	FAX	WEB
Germany	+49 (0) 91 26 / 25 87-25	www.physiomed.de

DEALER MARK / STAMP

TECHNOLOGY FOR THERAPY



