


# DEEP OSCILLATION®

DEEP OSCILLATION® is a unique international patented non-invasive, non-traumatic therapy method. Its special structure allows you to create biologically effective oscillations in the treated tissue using electrostatic attraction and friction. In contrast to other therapies, these pleasant oscillations have a gentle and deep-acting effect on all tissue components (skin, connective tissue, subcutaneous fat, muscles, blood and lymph vessels).

For the treatment, the patient holds a titanium neutral element loosely between the fingers. The pleasant therapy effect of DEEP OSCILLATION® is created beneath special gloves of the therapist or a specially designed hand applicator (second contact) circling over the tissue.

The exceptionally gentle action, early implementation and specific clinical effects make DEEP OSCILLATION® a unique therapy option, which is being used increasingly in medical fields that are not typically focused on by conventional physical therapies.



The following effects of DEEP OSCILLATION® are clinically proven:

- High pain-reducing potential (both for acute traumatic and chronic pain syndromes)
- Prevention and reduction of secondary and primary lymphoedema
- Muscular relaxation, movement-promoting effect, mobilisation
- Promotes wound healing, particularly for burns

ADRESS  
PHYSIOMED ELEKTROMEDIZIN AG  
Hutweide 10  
91220 Schnaittach  
Germany

PHONE  
+49 9126 2587 - 0  
  
FAX  
+49 9126 2587 - 25

E-MAIL  
info@physiomed.de  
  
WEB  
www.physiomed.de



# DEEP OSCILLATION® in Lymphology

2023-12. Errors and alterations reserved.



PHYSIOMED®  
technology for therapy



PHYSIOMED®  
technology for therapy

# Lymphoedema and lipoedema

The use of DEEP OSCILLATION® leads to volume reduction in primary and secondary lymphoedema. MLD and CPE can be intensified by DEEP OSCILLATION®. It is also increasingly and successfully used in the treatment of lipoedema.

DEEP OSCILLATION® has a loosening/reducing effect on fibrotic involvement. In particular, incipient fibrosis can be treated effectively. Various parameters concerning the quality of the skin are significantly improved. In addition, DEEP OSCILLATION® has the effect of pain reduction.



# Pre- and postoperative therapy

With DEEP OSCILLATION®, oedema can be eliminated much faster pre- and postoperatively than with conventional therapies. Due to extremely early application, pain is sustainably dampened and mobilisation is promoted. The quality of the tissue improves.

DEEP OSCILLATION® is therefore often routinely used as an adjuvant pre- and postoperative therapy, for example in liposuction. Because of its extremely gentle therapeutic effect, DEEP OSCILLATION® is used as a „first-line“ therapy after mastectomy, plastic surgery, osteosynthesis operations, endoprosthesis operations and caesarean section.



The clinical effects of DEEP OSCILLATION® are documented in over 60 studies. Here you can find more information and read about the individual studies.



Simply scan the QR code with your smartphone!

# DEEP OSCILLATION® Evident

Multifunctional intensity controls allowing for fast intensity reduction and quick switching between channels

Treatment index with filtering function



DEEP OSCILLATION® Evident allows you to begin treatment with maximum speed and ease: direct selection of deep oscillation parameters, over the indication index or program memory. During treatment you have a constant overview of all values and timers. With decades of proven use, the PHYSIOMED one-button operation allows fast intuitive operation in combination with the touch screen.

Deep oscillation characteristics can be adjusted precisely to the desired treatment over the individual setting of frequency (or frequency band), duty cycle, burst function and therapy duration.

The extensive indication list with practical information, graphics, animations and patient data-base with „potpourri“ function additionally ease use.

Up to three indication menus can be selected for DEEP OSCILLATION® Evident:

- CLINICS (for use in clinics, hospitals and physiotherapy practices)
- SPORTS (for use in sport medicine centres and clinics)
- AESTHETICS (for use in aesthetic medicine)

# DEEP OSCILLATION® Personal

Generous colour display for a good overview

Intuitive PHYSIOMED one-button operation

Treatment index for quick location of the desired treatment proposal featuring informative therapy notes and graphics for treatment



DEEP OSCILLATION® Personal can be used to read and retrieve program sequences previously written on the DEEP OSCILLATION® therapy card with the DEEP OSCILLATION® Evident unit in the clinic. This enables personalised mobile therapy on the ward and continuation of treatment in the patient's home according to the desired specifications.

This way you can guarantee the sustainability of your treatments until the next appointment in your clinic or treatment room.

DEEP OSCILLATION® Personal is also available with the Basic and Pro therapy cards. The Basic version offers a wide selection of the most important indication settings with pre-programmed treatment parameters, which usually consist of different treatment sections of different frequencies.

In addition, the Pro therapy card allows all parameters to be set manually as required for up to three therapy sequences of a treatment.